

Technical data sheet

Collection dummy for boars

Product description

Height adjustable dummy by means of a gas pressure spring. Heavy steel body coated with robust and removable PUR plastic cover to provide hygienic surface, removable mounting element, foot rests and base with locking pin. A later attachment of BoarMatic mechanics is possible.

Application

Designed for hygienic and efficient semen collection with maximum comfort for the boar.

Properties

Dimensions

Dummy / 11 inch	(l) 860 x (w) 285/500 x (h) 532-703 mm
Dummy / 9 inch	(l) 860 x (w) 230/450 x (h) 532-703 mm

Weight

Weight Dummy / 11 inch	80 kg
Weight Dummy / 9 inch	76 kg

Order information


Boar dummy with 11 inch top	REF. : 11100/0207
Boar dummy with 9 inch top	REF. : 11100/0107

Replacement parts

Mounting element	REF. : 5011100/0232
Plastic cover	REF. : 5011100/0233

Accessories and Consumables

Semen collection cup complete	REF. : 11123/1000
Blue Bag	REF. : 11123/0100
Semen collection bag 1 l	REF. : 11123/3010
Semen filter	REF. : 11123/3020
Rubber mat for boar dummy	REF. : 11100/0040
Mobile support for boar dummy	REF. : 11100/0110



Mounting element



Plastic cover



Mobile support



minitube