



# SEMEN COLLECTION- HANNOVER PHANTOM



The Hannover phantom has been designed to combine highest security for the stallion and the collection technician with the best possible comfort.

## Your benefits

- Adjustable in height and inclination to meet needs of different breeds and preferences.
- Covered with reinforced synthetic coating for a very high level of hygiene: the phantom can easily be washed and effectively be disinfected.
- Well padded to ensure highest comfort even to heavily stimulated stallions.
- Contoured body style makes it possible to hold and offer the AV at a most natural position.
- Easy to install and very stable: two tubular legs fixed on the floor with four screws each ensure highest security for operators and horses.
- Delivered with a removeable bite belt.

The Hannover phantom is available in two models: with and without hydraulic height adjustment for both legs. The hydraulic height adjustment ensures low fatigue operation. Both feet can be lifted or lowered about 40 cm either individually or simultaneously by means of a hand pump.

## Size of the phantom fits for all breeds:

Body: 200 x 60 x 55 cm (length x width x height)  
Height of the dummy: 1.30 m to 1.70 m approx.

## Hannover phantom

<u>with hydraulic height adjustment</u>	<b>REF. : 11200/0000</b>
<u>without hydraulic height adjustment</u>	<b>REF. : 11200/1000</b>



## Accessories

Protection covers: collection side closed with belts for attachment

**Leather blanket** (without illustration) [REF. : 11200/2501](#)

**Cover, synthetic material** [REF. : 11200/2601](#)



Bite belt, exchangeable, with adjustable fixing strips

**Bite belt** [REF. : 11200/2100](#)



Stretch foil for protection of the phantom collection end, length 300 m:

- For each collection a layer of the stretch foil is placed over the dummy mare to fulfil highest hygiene requirements.
- Easy application by wrapping it around the dummy, no additional fixing necessary.
- Less contamination of the dummy, cleaning is facilitated.
- After finishing all collections of the day all foil layers are removed in total.

**Stretch foil** [REF. : 32000/0010](#)

