



EASYCODER




The innovation in semen laboratories: printing of straws without ink

Your benefits

- Unmatched neat and clear printing on 0.5 and 0.25 ml straws
- Up to 3 lines with a total of about 200 characters, even on Ministraws
- Printing of barcode signs, logos and specific characters
- Immediately non-abrasive print
- Control via PC by a clear and user-friendly Minitube software
- Sequential numbers can be printed on straws
- Adjustable character fonts and graduations from micro to macro
- Freely eligible number of lines and of printing positions
- Clean and odourless working without ink and without solvents
- Easy handling, little maintenance required
- Easy and quick conversion between different straw sizes

EasyCoder printer

REF. : [13038/0000](#)



Technical features

- Thermal transfer printer with automatic straw supply
- Output: 3 600 straws per hour
- Windows based Minitüb software with data base for storing names and numbers
- Data transfer by serial port or USB adapter (included)
- Individual composition of printing lines possible, different configurations
- Dimensions (l x w x h): 690 x 370 x 310 mm (PC not included)
- Power supply: 110-230 V AC, 50/60 Hz

Accessories

Color ribbon for straw printer EasyCoder 2.0,
200m/roll, black

REF. : [13038/0010](#)