



PRINTERS AHDR AND HDR



AHDR – automatic printer

- For 0.5 ml or 0.25 ml straws and for Minitubes
- With counter and integrated motor
- Enables continuous printing of straws without manual effort
- Capacity: approx. 5 000 straws/hour

AHDR


REF. : [13142/1000](#)

HDR – manual printer

- The hand wheel operated HDR is a fast, clean and simple device for printing on straws (0.5 ml or 0.25 ml), as well as on Minitubes
- The device can easily be adjusted from one type of straw to the other
- Capacity: approx. 3 000 straws/hour

HDR

REF. : [13140/0000](#)



Accessories

Set of alpha-numerical characters for manufacturing of printing stamps

for medium straws 0.5 ml and Minitubes [REF. : 13052/0000](#)

for mini-straws 0.25 ml [REF. : 13052/0010](#)

Conversion kit for HDR/AHDR

from 0.5 ml to 0.25 ml straws [REF. : 13140/0525](#)

Printing ink ER, 40 ml, black [REF. : 13056/0033](#)

Thinner for ER printing ink, 1000 ml [REF. : 13056/9000](#)

Doublestick tape, 0.12 mm [REF. : 13050/0010](#)

