



MINITUBE SEMEN STRAWS



For over 25 years, the Minitube 0.25 ml and 0.5 ml semen straws have been used worldwide as a standard for the freezing of semen from bulls, stallions, boars, and many other species.

Minitube's quality control process, production techniques, and material selection ensure for a reliable product that consistently exceeds industry standards.

Minitube straws are available in a variety of transparent and full colors.

EcoStraws feature a plug with a shorter inner wick section that absorbs 3 - 4 % fewer sperm cells than the corresponding symmetrical plugs of standard semen straws. They are ideal for overnight equilibration of semen at refrigerated temperatures.

EcoStraws are used in exactly the same way as standard semen straws:

- identical processes for filling,
- identical processes for sealing,
- identical processes for printing,
- identical processes for insemination.



EcoStraw benefits

- Better utilization of available sperm cells per ejaculate because fewer sperm cells are lost in the plug
- Yield a larger number of semen straws per ejaculate
- Better use of valuable genetics without adding additional costs
- Ideal for overnight equilibration and fresh semen preservation



Colour	Mini straw 0.25 ml	EcoStraw 0.25 ml	Fresh semen straw 0.25ml	Medium straw 0.5 ml	EcoStraw 0.5 ml
Transparent	REF.: 13407/0010	REF.: 13407/3010	REF.: 13407/2010	REF.: 13408/0010	REF.: 13408/3010
Transparent – blue-white wick	REF.: 13407/0011	REF.: ---	REF.: ---	REF.: ---	REF.: ---
Red	REF.: 13407/0040	REF.: 13407/3040	REF.: ---	REF.: ---	REF.: ---
Red – transparent	REF.: 13407/0044	REF.: 13407/3044	REF.: 13407/2044	REF.: 13408/0044	REF.: 13408/3044
Green	REF.: 13407/0050	REF.: 13407/3050	REF.: ---	REF.: ---	REF.: ---
Green – transparent	REF.: 13407/0054	REF.: 13407/3054	REF.: 13407/2054	REF.: 13408/0054	REF.: 13408/3054
Blue – transparent	REF.: 13407/0064	REF.: 13407/3064	REF.: 13407/2064	REF.: 13408/0064	REF.: 13408/3064
Grey	REF.: 13407/0070	REF.: 13407/3070	REF.: ---	REF.: ---	REF.: ---
Grey – transparent	REF.: 13407/0074	REF.: 13407/3074	REF.: ---	REF.: 13408/0074	REF.: 13408/3074
Purple	REF.: 13407/0080	REF.: 13407/3080	REF.: ---	REF.: ---	REF.: ---
Yellow	REF.: 13407/0090	REF.: 13407/3090	REF.: ---	REF.: ---	REF.: ---
Yellow – transparent	REF.: 13407/0094	REF.: 13407/3094	REF.: 13407/2094	REF.: 13408/0094	REF.: 13408/3094
Pink	REF.: 13407/0100	REF.: 13407/3100	REF.: ---	REF.: 13408/0100	REF.: 13408/3100
Turquoise	REF.: 13407/0130	REF.: 13407/3130	REF.: ---	REF.: 13408/0130	REF.: 13408/3130
Orange	REF.: 13407/0140	REF.: 13407/3140	REF.: ---	REF.: 13408/0140	REF.: 13408/3140
Orange – transparent	REF.: 13407/0144	REF.: 13407/3144	REF.: ---	REF.: 13408/0144	REF.: 13408/3144
Salmon	REF.: 13407/0160	REF.: 13407/3160	REF.: ---	REF.: ---	REF.: ---
Beige	REF.: 13407/0180	REF.: 13407/3180	REF.: ---	REF.: ---	REF.: ---
Pistachio green	REF.: 13407/0200	REF.: 13407/3200	REF.: ---	REF.: 13408/0200	REF.: 13408/3200
Pistachio green transparent	REF.: 13407/0204	REF.: 13407/3204	REF.: ---	REF.: 13408/0204	REF.: 13408/3204
White	REF.: ---	REF.: ---	REF.: ---	REF.: 13408/0210	REF.: 13408/3210