



# BOVINE ARTIFICIAL INSEMINATION



## QuickLock AI gun

**The AI gun for AI success.**

- All-purpose gun for all types of straws and sheaths, split or non-split
- Stainless steel and sterilisable parts only
- One hand only sheath fixation and release
- Unique design for safe and centered straw alignment every time – no waste of straws
- Separate straw release for control of information after AI
- No loose parts – safety catch prevents the plunger from sliding out
- Identification with different coloured stylets
- All parts can be replaced

## QuickLock Classic

Sheath fixation and release very easy by just pressing a knob

**QuickLock Classic** REF. : 17025/0000

## QuickLock 2000

- Ergonomic design
- Only 1 cm of the sheath length is required for fixation

**QuickLock 2000** REF. : 17026/0000

## QuickLock ERGO

- Ideal for older cows due to increased gun length
- Smallest ergonomic handle available
- Easiest gun as far as cleaning is concerned

**QuickLock ERGO** REF. : 17024/0000

## QuickLock stylets

Available in different colours.


## Universal sheaths

- Suitable for all types of semen straws
- Absolutely tight fit with all types of straws
- Gentle insemination due to a soft shape of the tip
- 50/bag, 1000/box
- Transparent adapters for better visual control of straw position

**Universal sheath split** (transparent-green) REF. : 17007/0001

**Universal sheath non-split** (transparent-blue) REF. : 17007/0000

**Hygienic sheath**, 80/roll REF. : 17028/0100





## Thawing unit MT 30/54

- Electronic temperature control
- Central, built-in thawing chamber with lift
- Accuracy: 0.2°C
- Power supply: 12 V DC (car battery)

35°C	REF. : 17044/0035
37°C	REF. : 17044/0037
38°C	REF. : 17044/0038
50°C	REF. : 17044/0050



Transformer for thawing unit 12 V/220 V REF. : 17165/1395

Thermometer for thawing unit REF. : 17049/0000

## Insemination kit

Insemination kit complete: heavy duty case, electronic thawing unit MT 30/54, QuickLock, 80 Universal sheaths, ReproJelly lubricant, speculum Varikon, 1 box AI gloves, MiniCutter, tweezers and clamp

Insemination kit REF. : 17050/0060

Components also available individually:

Insemination case	REF. : 17050/0001
Insemination gloves, 100/box	REF. : 17080/0015
Varikon speculum	REF. : 17064/0380
Straw tweezers, angular shape	REF. : 17060/0000
ReproJelly, 100 ml	REF. : 11907/0100
MiniCutter	REF. : 17062/0010
Clamp	REF. : 17063/0000

Not included in kit:

Transformer for thawing unit REF. : 17165/1395

## QuickLock Heater 2.0

battery operated, with carrying strap

- AI guns are transported at a constant temperature
- Thawed semen will not lose temperature in pre-warmed QuickLock AI gun
- Very robust and easy to clean material, bag can be washed in a washing machine
- Strap allows different ways of carrying the bag
- Rechargeable battery, works for approx. 8 hours
- Temperature setting : +35°C, other values optional
- Free choice of charger, either 220 V or 12 V/ 24 V



QuickLock Heater 2.0

REF. : 17028/0210

## Accessories

Charger for recharging of battery set:

Charger for wall socket 100-240 VAC, 50/60 Hz REF. : 17028/0220

Charger for car, 12 V/24 V REF. : 17028/0221

Disposable plastic bag

to protect the inside, 100/bag REF. : 17028/0201

Additional rechargeable battery set

REF. : 17028/0219

