



# QUICKLOCK HEATER 2.0




## Your benefits

- AI guns are transported at a constant temperature
- Temperature of thawed semen will not drop in pre-warmed QuickLock AI gun
- Very robust material, easy to clean, bag can be washed in a washing machine
- Strap allows different ways of carrying the bag
- Rechargeable battery, works for approx. 8 hours, dependent on ambient temperature
- Temperature setting range: +35°C, other values upon demand

## QuickLock Heater 2.0

with battery-set and disposable plastic bag [REF. : 17028/0210](#)



Insemination guns  
not included



## Accessories

### Charger for recharging of battery set:

Charger for wall socket 100-240 VAC, 50/60 Hz [REF. : 17028/0220](#)

Charger for car, 12 V/24 V [REF. : 17028/0221](#)

### Option:

Additional rechargeable battery set [REF. : 17028/0219](#)

QuickLock heater 2.0, set to 25°C for ET [REF. : 17028/0225](#)

### Disposable plastic bag

to protect the inside, 100/bag [REF. : 17028/0201](#)



minitube

[www.minitube.com](http://www.minitube.com)