



# CLEARGLIDE®



## HIGH QUALITY INSEMINATION CATHETER WITH UNIQUE POLY GEL TIP

The ClearGlide® tip locks firmly into the sow's cervix due to the unique flexible collar that flips back during insertion but once in place, it naturally creates a firm grip. The Poly Gel tip conforms naturally to the shape of the cervix. The surface of the nonabsorbent material has a special powder coating that creates a slippery film when it comes into contact with moisture, so no extra lubricant is necessary.

Due to the unique tapered shape and sealing flange of the tip, the ClearGlide® offers superior locking of the catheter into the cervix and therefore reduces semen back-flow down to a minimum. Recessed shaft and superior gliding capacity make insertion and removal effortless and improve the animal well-being.




### Your Benefits

- ClearGlide® tip locks firmly into the sow's cervix due to the unique flexible collar that flips back during insertion.
- Superior gliding capacity and the tapered shape allow gentle insertion and removal.
- Flexible Poly Gel tip adjusts perfectly to the cervix.
- Semen back-flow is reduced down to a minimum.
- The recessed shaft prevents irritation and injury to the cervix.
- 'All in one' moulded manufacture provides firm fixation of tip to shaft.
- The high standard of automatic manufacture guarantees quality and biosecurity.





## Superior lock of ClearGlide® in the cervix



### Product features

- New material with unique molecular structure
- Poly Gel tip has different hardness zones – soft in the front – more dense in the rear
- Tip is moulded to the recessed shaft
- Material offers extraordinary gliding capacity
- The highest hygiene standard available in the SafeBlue packaging

### SafeBlue ClearGlide®

ClearGlide® is available in the hygienic SafeBlue presentation: top quality insemination catheters, single packed and prelubed. Better insemination hygiene improves reproduction results.



**ClearGlide® with handle, 25/bag,  
500/box**

REF. : [17106/4026](#)

**ClearGlide® without handle, 25/bag,  
500/box**

REF. : [17106/3026](#)

**SafeBlue ClearGlide®  
with handle, 100/package**

REF. : [17106/4005](#)

**SafeBlue ClearGlide®  
without handle, 100/package**

REF. : [17106/3005](#)