



AI TRAINING MODEL "SUZY"



A new approach to artificial insemination training

"Suzy" is an artificial sow for training AI technicians in a safe and convenient environment, before they start working with real animals.

The training model consists of a plastic body and silicone internals. It provides a **most realistic gynaecological tract** (GT) that can be removed to demonstrate anatomy and the AI process as an initial step of the training session.

"Suzy's" GT features a cervix, similar to the natural one, for practising the lock of the AI catheter in the cervix. A sound will confirm success when the introduced catheter reaches the correct position in the cervix („oink"). A replica of the urethra opening trains the correct introduction of the AI catheter avoiding positioning it in the urethra. Inseminating "Suzy" with a „real" semen dose is possible, as the GT is provided with a drain. Both, **traditional** and **post-cervical AI** can be practised when using the PCAI genital tract. This tract features a clear tube, where the trainer can overview the advance of the inner cannula of the PCAI catheter.

The artificial sow provides flexibility around timing and location of AI training courses, and also lowers the risk of injury to animals by unexperienced trainees. Training units can be repeated as often as necessary. All these benefits help to provide **optimal training** and ultimately an improved performance of applicants of AI in sows.

Product features

- Anatomically correct sow model including full reproductive tract and cervix which provides an authentic feeling during AI training
- Modular approach providing different training levels for quick learning success
- Passage of inner cannula of PCAI catheter is very natural through noticeable cervix cushions
- Correct introduction of AI catheter and inner cannula is confirmed by oinking sound
- „Real", complete insemination process can be performed
- Compact and light-weight training model for easy transport and cleaning



Increase the efficiency of your AI training courses and contribute to the enhancement of animal welfare!



Order information

Artificial sow "Suzy" - with complete gynaecological tract for traditional cervical insemination

REF. : [22400/0100](#)

Artificial sow "Suzy" - with 2 complete gynaecological tracts for traditional and post-cervical insemination

REF. : [22400/0101](#)

Starter-kit includes

Ref. 22400/0100:

100 AI catheters and 2 tubes of ReproJelly

Ref. 22400/0101:

100 AI catheters, 25 PCAI catheters and 2 tubes of ReproJelly



Spare parts

Outer part with rectum, vulva and tail

REF. : [22400/0105](#)

Genital tract for traditional and post-cervical insemination with Oinkmaster, vaginal tube and cervix

REF. : [22400/0106](#)

