



# MINITUBE ULTRASOUND SCANNER



## Portable ultrasound scanner for transabdominal pregnancy detection in sows and small ruminants

### Your benefits

- compact and light weight
- portable, battery operated
- robust, electric probe
- very easy to clean
- integrated auto-measure function to easily measure back fat thickness

Stable carrying case, electronic 3.5 MHz convex sensor (70° forward), freeze-function, carrying strap, two rechargeable batteries, battery charger, USB cable to transmit photos

**Ultrasound scanner Minitube** REF.: **23500/1000**



### Accessories

**Spare battery set (2 pc.)** REF.: **5023500/1010**

**ReproJelly, lubricant, 100 ml** REF.: **11907/0100**

### Technical details

Monitor size:	5.0" TFT-LCD
Display modes:	B, B+B
Image gray scale:	256 levels
Measure functions:	Distance, circumference, area, volume
Detection depth:	up to 240mm
Battery capacity:	3000mAh, 7.4V
Power consumption:	8W
Weight (incl. probe):	854 g



Pregnant sow, 22 days



Pregnant sow, 29 days



Pregnant sow, 40 days



Pig foetus, 60 days