



MINITUBE ULTRASOUND SCANNER 2.0



PORTABLE ULTRASOUND SCANNER FOR PREGNANCY DETECTION

Your benefits

- Good image quality at an affordable price
- Easy documentation of the analyses through labelling of images and recording of short video sequences
- A variety of probes can be employed (see accessories)
- Compact, light weight and very robust
- You can use the equipment while the extra batteries are recharging

Comes with two rechargeable batteries, charger, carrier bag, contact gel and heavy duty carrying case.

Ultrasound scanner Minitube 2.0 REF.: **23500/1005**

Accessories

6.5 MHz Linear Rectal Probe 1

Range: 6.5 / 7.5 / 8.5 MH REF.: **23500/1106**

3.5 MHz Convex Probe 2

Range: 2.5 / 3.5 / 5.0 MH REF.: **23500/1105**

Technical details

Monitor size:	7.0" TFT-LCD
Scanner dimensions:	216 x 147 x 27 mm
Display modes:	B, B+B, B+M, M, 4B
Image gray scale:	256 levels
Image multiplying factor:	×1.0, ×1.2, ×1.5, ×2.0
Measure functions:	Distance, circumference, area
Port:	USB 2.0
Battery capacity:	3000 mAh
Power consumption:	7 W
Weight (excl. probe):	900 g
Voltage:	100 V-240 V
Automatic probe identification	

www.minitube.com

